


## Exercise 10-4: Programmatically Create Filenames Based on the Current Timestamp

### Goal

- Explore example code that programmatically creates filenames based on the current date and time.

### Examine Code to Generate a Timestamped Filename in Relative Path

- Open the following project: C:\Exercises\LabVIEW Core 1\ Timestamped Filepath\Timestamped Filenames.lvproj.
- From the **Project Explorer** window, open the Generate Timestamped Filepath VI.
- Examine the block diagram.



- Format Date/Time String function** – Converts a timestamp of the current time or numeric value into a string that displays the corresponding time. Using the time format string in this block diagram, this function will return a string similar to the following:

<year><month><day>\_<hour><minute><second>\_Acquired Data.txt

Example: 201901014\_101530\_Acquired Data.txt

- Application Directory VI** – Returns the path to the directory containing the current project (.lvproj).
- Build Path function** – Creates a new path by appending a filename or relative path to an existing path. In this VI, the output of this function will be a timestamped filename in the directory containing the current project. (For example, C:\Exercises\LabVIEW Core 1\Timestamped Filepath\ 20170601\_101530\_Acquired Data.txt)




**Note:** Refer to the *LabVIEW Help* for how to use the **time format string** input

4. Test the code.
  - Run the VI multiple times.
  - On the front panel, notice that the generated filename changes to match the current date and time.
  - Notice that the filepath directory is the same as the directory containing the current project file (.lvproj).

### Examine Timestamped Filename Code in a Data Logging VI

1. From the **Project Explorer** window, open the [Timestamped File] High-Level Write to Text File VI.
2. Examine the block diagram. Notice how the timestamped filepath code passes the timestamped filepath to the **filepath input** of the file Write Delimited Spreadsheet VI.



3. Test th VI.
  - Run the VI.
  - Navigate to the log file in the C:\Exercises\LabVIEW Core 1\Timestamped Filepath directory. Notice that the filename indicates the date/time when the VI created that file.
  - Run the VI a couple more times.
  - Examine the additional log files that the VI generated in the C:\Exercises\LabVIEW Core 1\Timestamped Filepath directory. Notice that their filenames also indicate the date/time when the VI created those files.

## Your Turn

Modify the VI so that the VI creates a timestamped filename that looks similar to the following format and saves the file in the C:\Exercises\LabVIEW Core 1\ Timestamped File directory.

<year><month><day>\_<hour><minute><second>\_<your own custom filename>

Example: 20190511\_090030\_Batch ABC.txt

## On the Job

Would any of your applications benefit from creating a new timestamped filename every time you run your VI?

---

---

---

---

---

**End of Exercise 10-4**