

Exercise 2-2: Validating I/O (NI DAQ Device)

Goal

To verify that the inputs and outputs of your DAQ device behaves correctly.

Instructions


Examine the configuration of DAQ device:

1. Launch **NI MAX** from the Windows **Start** menu.
2. Expand **Devices and Interfaces** and select the **PCI-6221** or **USB-6212**.



Note that from here till the end of the course choose the NI-DAQ device that you have connected to your PC, e.g. if you are using USB-6212, then select it each time you see **PCI-6221** is mentioned.

3. Right-click the **device** and select **Test Panels** from the shortcut menu.
4. Validate the analog input of your device.
 - Make sure that the **Analog Input** tab is selected in the Test Panels dialog box.
 - Select **PCI-6221/ai2** from the **Channel Name** control, which corresponds to the ai2 pin of the device.
 - Change **Mode** to **Continuous** enables the measurement-specific settings, e.g. **Rate (Hz)**, and **Samples to Read**.
 - Click the **Start** button to start the acquisition.
 - Explore data on the chart, then click **Stop**.
 - Change **Samples To Read** to other values, e.g. 5000, 10000, etc.
 - Start the acquisition, and explore how it affects the acquired signal.



- Click the **Stop** button.
 - You can disable the **Auto-scale chart** option, and change the scaling of the x-scale to adjust the graph by double-clicking the **maximum** and **minimum** values of the scale and typing the appropriate ones.
5. Validate the analog output of your device.
 - Switch to the **Analog Output** tab.
 - Set the Channel Name control, e.g. PCI-6221/ao0.
 - You can either change Mode to Voltage Sinewave Generation or leave Voltage DC.
 - Set **Output Value (V)** to the constant voltage you want to output.
 - Click the **Update** button to start the generation.
 - Click **Stop**.
 6. Validate the digital I/O of your device.
 - Switch to the **Digital I/O** tab.
 - Select the **port** from the **Port Name** control. Port is a collection of lines, and each line can output a high or low voltage.
 - Select **All Output** in the Select **Direction** section.
 - Configure the output voltage in the **Select State** section by changing the line state.
 - Click **Start**.
 - Select **All Low** and click **Stop**.
 - Close the **Test Panels** dialog box.
 - Close NI MAX.

End of Exercise 2-2